


How to use a fentanyl test strip

What do I need?


- ✓ **1 shot glass, half-filled with water**



You can use any small container. Little ketchup cups work well! Whatever container you use, only use about a half (1/2) a shot glass of water (15mL). If you are testing a stimulant (like meth, cocaine, or MDMA) use a full shot glass of water instead.




- ✓ **fentanyl test strip**



- ✓ **very small amount of powder to test**


about the size of a match head or half a grain of rice



about this much!

How do I use the strip?

- 1 Put the powder in the glass of water and swirl it until it is completely dissolved. (If you want to test something that isn't a powder, crush it first.)




- 2 Hold the solid blue end of the test strip and dip the other end in the water.

- 3 Hold the strip in the water for about ten seconds.

You will know it's ready because the water will be sucked up into the white part of the strip.

- 4 Take the strip out of the water and wait sixty (60) seconds.

- 5 Use a bright light to look at the strip and count how many red lines there are. Even very light red lines count!



TWO (2) red lines means... **NO**, there is no fentanyl!



ONE (1) red line means... **YES, there is fentanyl!**

No drug is completely safe. Make sure someone is with you who can give naloxone (Narcan) or call for help if you need it.

- 6 Make sure you check the strip within ten minutes of doing the test. Throw the strip away when you are done.

